

Staining Index and Spill Over

on

LSRFortessa

Determined using anti-CD8 in an
Experiment for Analysing Lymphocytes.

FACS Core Facility, Aarhus University
2016

Index

1. Index
2. Aim for the study, Materials and Methods
3. Antibodies, fluorochromes, lasers, emission wavelengths
4. SOP for staining: Antibody staining of surface antigens on cells using a 96 well plate
- 5. The "Compensation Matrix" and PMT voltages**
- 6. Stain Index**
- 7. Fluorochrome Suggestions for 4-8 color panel**

Extras

8. Compensation Beads : Fluorescence histograms
9. Lymphocytes: Gating on the fluorescence histograms
10. Stain Index calculations
11. Lasers and Filters at the LSRFortessa (May 2016)
(See the newest version on www.facs.au.dk/Links/LSRFortessa)

Aim for this study

To get to know our specific flow cytometer – LSRFortessa – better, we decided to calculate Compensation Values and Stain Index for 25 different fluorochromes in an experiment optimized for running human peripheral blood lymphocytes.

It is important to realize, that the compensation values and stain index values are closely related to these settings.

Materials

Antibodies:

25 different fluorochromes conjugated to anti-human CD8 antibody, clone RPA-T8 except for BV605, PerCP-eFlour710, APC (clone SK1), and PE-Cy7 (clone HIT8a). (See page3 for details concerning antibodies, fluorochromes, lasers and emission wavelengths)

MNC separated from human peripheral blood.

Compensation Beads: eBioscience 01-1111-42

Other materials: See the procedure “Antibody staining of surface antigens on cells - using a 96 well plate”.

Methods

Lymphocytes and staining.

MNC were separated from human peripheral blood by centrifugation on Histopaque gradient (Sigma). Washed twice in PBS with 0.5% BSA. MNC were frozen in RPMI-1640 containing 10% DMSO and 20% FCS.

MNC were thawed, centrifuged and resuspended in PBS containing 0.5% BSA and 0.09% NaN₃ and stained as described in “Antibody staining of surface antigens on cells using a 96 well plate” (page4) One sample for each antibody-fluorochrome-conjugate and one unstained sample.

OneComp eBeads were stained, as the cells were, but no unstained sample was included, as the bead suspension contains negative as well as positive beads.

An experiment for running lymphocytes

was created by running unstained mononuclear cells on the LSRFortessa, gating the lymphocytes and setting the median fluorescence from each PMT above the lower limit for linearity. Trying to avoid negative values after compensation, voltages were set a bit higher to achieve medians around 100-200.

Compensation beads and MNCs were analyzed in this experiment

Compensation Values

were determined using the beads FCS files and FlowJo Software.

As we wanted to calculate compensation for up to 3 fluorochromes per detector, we had to set up 3 different compensation matrices in FlowJo. After calculating the 3 matrices the values were pooled in one excel spreadsheet.

Stain index

were calculated as $(\text{Median}_{\text{pos}} - \text{median}_{\text{neg}}) / (2 \times \text{rSD}_{\text{neg}})$ after analyzing the lymphocyte gate from the MNC FCS files in FlowJo software.

Results can be seen at pages 5 The Compensation Matrix and 6 Stain Index

Antibody overview

| Laser | Well | µl/test ab | FL | klon | Cat # | Lot # | Exp date | Comp. matrice nr. * | Filter | Filter | |
|-------|------|------------|-----|-----------------|--------|------------------|-------------|---------------------|--------|--------|-----|
| A | A1 | 5ul | CD8 | V450 | RPA-T8 | BD560347 | 3039752 | 2015-01 | 1 | 442/46 | V1 |
| | A2 | 5ul | CD8 | BV421 | RPA-T8 | BD562429(25test) | 4290912 | 2018-10 | 2 | 442/46 | V1 |
| | A3 | 5ul | CD8 | V500 | RPA-T8 | BD560775 | 3235542 | 2015-03 | 1 | 525/50 | V2 |
| | A4 | 5ul | CD8 | BV510 | RPA-T8 | BD563256 | 4080668 | 2015-07 | 2 | 525/50 | V2 |
| | A5 | 5ul | CD8 | eVolve605 | RPA-T8 | eBio83-0088-41 | E24358-102 | 2015-11 | 1 | 610/20 | V3 |
| | A6 | 5ul | CD8 | BV 605 | SK1 | BD564116 | - | - | 2 | 610/20 | V3 |
| | A7 | 5ul | CD8 | BV650 | RPA-T8 | BD563822 | 4115843 | 2017-09 | 1 | 660/20 | V4 |
| | A8 | 5ul | CD8 | BV711 | RPA-T8 | BD563676 | 4108877 | 2018-02 | 1 | 711/25 | V5 |
| | A9 | 5ul | CD8 | BV786 | RPA-T8 | BD563824 | 5098928 | 2016-07 | 1 | 785/80 | V6 |
| B | B2 | 5ul | CD8 | PerCP-eFlour710 | SK1 | eBio46-0087-41 | E10832-1633 | 2017-02 | 2 | 710/50 | B2 |
| | B1 | 5ul | CD8 | PerCP-Cy5.5 | RPA-T8 | BD560662 | 4339888 | 2016-10 | 1 | 710/50 | B2 |
| | B12 | 5ul | CD8 | BB515 | RPA-T8 | BD564527 | 4255561 | 2015-12 | 3 | 530/30 | B1 |
| | A11 | 5ul | CD8 | AF488 | RPA-T8 | BD557704 | 4171577 | 2015-10 | 2 | 530/30 | B1 |
| | A10 | 20ul | CD8 | FlTC | RPA-T8 | BD561947 | 4028695 | 2015-12 | 1 | 530/30 | B1 |
| | B3 | 20ul | CD8 | PE | RPA-T8 | BD561949 | 43044878 | 2020-08 | 1 | 586/15 | GG1 |
| C | B4 | 5ul | CD8 | PE-eFlour610 | RPA-T8 | eBio61-0088-41 | E19357-101 | 2016-02 | 1 | 615/20 | GG2 |
| | B5 | 5ul | CD8 | PE-CF594 | RPA-T8 | BD562311 | 3346941 | 2015-11 | 2 | 615/20 | GG2 |
| | B6 | 10ul | CD8 | PE-Cy5 | HIT8a | BD555636 | 2097558 | 2015-03 | 1 | 660/20 | GG3 |
| | B7 | 5ul | CD8 | PE-Cy7 | RPA-T8 | BD560917 | 3287705 | 2015-04 | 1 | 780/60 | GG5 |
| | B8 | 10ul | CD8 | APC | RPA-T8 | BD555369 | 12780 | 2014-03 | 1 | 670/30 | R1 |
| D | B9 | 5ul | CD8 | AF647 | RPA-T8 | BD557708 | 3252624 | 2015-04 | 2 | 670/30 | R1 |
| | B10 | 0.2ul | CD8 | AF700 | RPA-T8 | BD557945 | 3316609 | 2015-10 | 1 | 730/45 | R2 |
| | B11 | 5ul | CD8 | APC-eFlour780 | RPA-T8 | eBio47-0088-41 | E08444-1634 | 2016-09 | 1 | 780/60 | R3 |
| | B12 | 5ul | CD8 | APC-H7 | SK1 | BD560273 | 4238532 | 2016-02 | 2 | 780/60 | R3 |
| E | C1 | 5ul | CD8 | APC-Cy7 | RPA-T8 | BD557760 | 4045963 | 2015-12 | 3 | 780/60 | R3 |

* Which color goes to which Compensation Matrix calculation.

To calculate compensation matrices, we have to make 3 matrices as we test up to 3 colors in the same channel.

Antibody staining of surface antigens on cells using a 96 well plate

Notes

This is a general protocol, which does not concern specific antibodies, however it takes into consideration: Incubation, washing and fixation.

Consider blocking before staining.

In this protocol the incubation is performed at RT. If you want to incubate on ice the incubation time should be expanded to 30-60 min and the centrifuge should be cold.

In each well there will be approximately 15µl left after pouring off your supernatant.

MATERIAL

1. Stain buffer: PBS pH 7.4 with 0.5% BSA and 0.09% Na-azide
2. Fluorochrome conjugated antibodies and maybe isotype controls
3. Cells at $1-5 \times 10^6$ /ml in stain buffer
4. A round bottom microwell plate (not for cell culturing)
5. Fixation buffer: PBS pH 7.4 with 0.9% formaldehyde

PROCEDURE

Be cautious at any time not to mix fluid from one well to another when pouring fluid off. Even a tiny bit of antibody can stain your cells.

1. Make a scheme showing which antibodies in which wells. If viability test is not a part of your panel you should make a sample just with your viability marker
2. Adjust your cell suspension to $1-5 \times 10^6$ /ml in stain buffer. Blocking could be performed now
3. Add antibodies to the wells
4. Transfer 100µl cell suspension to each well and pipette up and down 5 times to mix cells and antibodies.
5. Incubate in the dark 15-30 min at RT. Some cell types will require incubation 30-60 min at 4 C⁰
6. Add 100 µl washing buffer to each well.
7. Centrifuge the plate at 350g for 2 min at RT (or cold)
8. Place a thick tissue flat on the table
9. Pour off the supernatant in one sliding movement and press the plate briefly against the tissue before you turn the plate bottom downwards again
10. Loosen the cell pellet by a gentle knocking on the side of the plate
11. Add 200µl stain buffer to each well using a multichannel pipette, resuspend cells using the pipette.
12. Repeat from 7-10 twice
13. Fix and resuspend the cells by adding 200µl of fixation buffer to each well and pipette up and down 5 times
14. Place a lid on the plate and keep it cold and dark until analysis
15. The use of tandem dyes will limit the storage time to 24 hours. We recommend to analyse as quickly as possible.

LSRFortessa Compensation Matrix for "Lymphocyte Settings" See these settings below **More than one stain per PMT is shown!**
 Staining performed May 2015 and Feb 2016

Spill over from row to column


| PMT | V1 | V2 | V3 | V4 | V5 | V6 | B1 | B2 | YG1 | YG2 | YG3 | YG5 | R1 | R2 | R3 | % comp. |
|------------------|----------|----------|----------|----------|----------|----------|-------------|-------------|-----------|-----------|-----------|-----------|------------|------------|------------|-----------|
| Band Pass filter | V 442_46 | V 525_50 | V 605_15 | V 660_20 | V 711_25 | V 785_60 | Blue 530_30 | Blue 710_50 | YG 575_15 | YG 615_20 | YG 660_20 | YG 780_60 | Red 670_30 | Red 730_45 | Red 780_60 | color key |
| V450 | | 25,6% | 4,0% | 0,8% | 0,3% | 0,1% | 0,0% | 0,0% | 0,0% | 0,0% | -0,1% | 0,0% | -0,3% | -0,1% | 0,0% | 0-5 |
| BV421 | | 11,3% | 1,6% | 0,4% | 0,1% | 0,1% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 5-10 |
| V500 | 3,3% | | 53,7% | 14,7% | 4,9% | 1,8% | 5,9% | 0,3% | 0,0% | 0,0% | 0,0% | 0,0% | -0,2% | -0,1% | 0,0% | 10-20 |
| BV510 | 7,5% | | 99,0% | 36,3% | 15,2% | 7,3% | 0,7% | 0,1% | 0,0% | 0,0% | 0,0% | 0,0% | -0,1% | 0,0% | 0,0% | 20-30 |
| eVolve605 | 0,0% | 0,0% | | 1,4% | 0,1% | 0,1% | 0,0% | 0,0% | 1,5% | 5,0% | 0,2% | 0,0% | -0,1% | 0,0% | 0,0% | 30-40 |
| BV605 | 1,5% | 0,2% | | 49,4% | 18,0% | 7,8% | 0,0% | 1,2% | 5,1% | 8,5% | 12,4% | 1,2% | 0,0% | 0,1% | 0,0% | 40-60 |
| BV650 | 1,7% | 0,3% | 17,4% | | 40,1% | 13,5% | 0,0% | 0,8% | 0,1% | 0,9% | 7,1% | 0,6% | 8,4% | 8,3% | 1,9% | 60-80 |
| BV711 | 4,4% | 0,7% | 0,5% | 2,6% | | 53,5% | 0,0% | 4,2% | 0,0% | 0,0% | 0,1% | 0,4% | 1,0% | 25,4% | 7,7% | 80-100 |
| BV786 | 3,9% | 0,8% | 0,7% | 1,1% | 2,2% | | 0,0% | 0,1% | 0,0% | 0,0% | 0,0% | 0,5% | 0,2% | 2,0% | 7,1% | |
| FITC | 0,0% | 3,8% | 1,5% | 0,4% | 0,2% | 0,1% | | 1,8% | 0,0% | 0,0% | 0,0% | 0,0% | -0,1% | 0,0% | 0,0% | |
| AF488 | -0,6% | 0,9% | 0,3% | 0,1% | 0,0% | 0,0% | | 1,1% | 0,0% | 0,0% | 0,0% | 0,0% | -0,1% | 0,0% | 0,0% | |
| BB515 | -0,2% | 0,9% | 0,2% | 0,0% | 0,0% | 0,0% | | 0,9% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | 0,0% | |
| PerCP-Cy5.5 | 0,0% | 0,0% | 0,0% | 13,8% | 83,2% | 39,2% | 0,0% | | 0,0% | 0,0% | 8,2% | 11,5% | 17,7% | 74,5% | 23,7% | |
| PerCP-eFlour710 | 0,0% | 0,0% | 0,0% | 2,0% | 77,6% | 35,7% | 0,0% | | 0,0% | 0,0% | 1,1% | 9,2% | 2,9% | 93,8% | 24,7% | |
| PE | 0,0% | 0,1% | 12,8% | 3,2% | 1,1% | 0,3% | 0,9% | 9,9% | | 46,6% | 31,0% | 2,0% | 0,0% | 0,0% | 0,0% | |
| PE-eFlour610 | 0,0% | 0,0% | 23,7% | 6,6% | 1,8% | 0,6% | 0,1% | 15,7% | 40,3% | | 81,9% | 4,6% | 0,1% | 0,1% | 0,0% | |
| PE-CF594 | 0,0% | 0,0% | 24,1% | 7,8% | 3,6% | 1,5% | 0,1% | 33,2% | 19,1% | | 87,1% | 10,1% | 0,3% | 0,4% | 0,1% | |
| PE-Cy5 | 0,0% | 0,0% | 0,2% | 8,3% | 4,0% | 1,7% | 0,0% | 40,1% | 1,5% | 0,8% | | 13,0% | 33,7% | 31,9% | 8,6% | |
| PE-Cy7 | 0,0% | 0,0% | 0,1% | 0,0% | 0,1% | 17,0% | 0,1% | 1,0% | 0,7% | 0,4% | 0,3% | | 0,0% | 1,3% | 16,9% | |
| APC | 0,0% | 0,0% | 0,2% | 11,6% | 3,1% | 1,4% | 0,0% | 0,7% | 0,0% | 0,6% | 44,7% | 3,3% | | 82,9% | 20,8% | |
| AF647 | 0,0% | 0,0% | 0,0% | 0,3% | 0,1% | 0,1% | 0,0% | 0,3% | 0,0% | 0,0% | 8,0% | 1,2% | | 97,8% | 26,8% | |
| AF700 | 0,0% | 0,0% | 0,0% | 0,0% | 4,0% | 2,2% | 0,0% | 1,1% | 0,0% | 0,0% | 0,1% | 2,1% | 1,1% | | 25,9% | |
| APC-eFlour780 | 0,0% | 0,0% | 0,0% | 0,6% | 0,2% | 8,3% | 0,0% | 0,1% | 0,0% | 0,0% | 2,3% | 10,8% | 6,7% | 11,1% | | |
| APC-H7 | 0,0% | 0,0% | 0,0% | 0,1% | 0,0% | 8,6% | 0,0% | 0,0% | 0,0% | 0,0% | 0,2% | 9,8% | 0,7% | 5,0% | | |
| APC-Cy7 | 0,0% | 0,0% | 0,0% | 0,4% | 0,1% | 7,8% | 0,0% | 0,0% | 0,0% | 0,0% | 1,4% | 10,6% | 3,9% | 11,5% | | |

"Lymphocyte settings":

| Short name | V1 | V2 | V3 | V4 | V5 | V6 | B1 | B2 | YG1 | YG2 | YG3 | YG5 | R1 | R2 | R3 |
|------------------|----------|----------|----------|----------|----------|----------|-------------|-------------|-----------|-----------|-----------|-----------|------------|------------|------------|
| Band Pass filter | V 442_46 | V 525_50 | V 605_15 | V 660_20 | V 711_25 | V 785_60 | Blue 530_30 | Blue 710_50 | YG 575_15 | YG 615_20 | YG 660_20 | YG 780_60 | Red 670_30 | Red 730_45 | Red 780_60 |
| PMT Voltage | 460 | 480 | 620 | 630 | 680 | 710 | 500 | 600 | 520 | 540 | 660 | 660 | 700 | 720 | 720 |

Short names for the PMTs: The PMTs are named from 1 and forthgoing for each laser, where V is violet laser and YG is yellow-green laser
 The PMT YG4 (710/40) has not been used in this experiment

Stain Index, Fortessa, CD8, 20.05.2015



Fluorochrome/PMT suggestions for 4-8 color panel

Based on our overview "Lasers and Filters at the LSRFortessa" (page 11) and a glance to the Compensation Matrix from this experiment, we will set up suggestions for choosing fluorochromes/PMTs for multicolor panels.

These suggestions are NOT the only way to do it!!

You should still remember to take into account the abundance of the target antigen, when you choose fluorochrome. Choose a fluorochrome with a high Stain Index for a weakly expressed antigen and vice versa.

An example: If you look for a fluorochrome for V1, you could choose BV421 for a weakly expressed antigen and V450 for a highly expressed antigen.

You can use more colors, if you have subpopulations, which do not express all antigens in your panel.

In stead of proposing fluorochromes for the panels, we will propose PMTs to use. (see page 3 or 5 to understand the abbreviations)

PMT suggestions

4 colors: V1, B1, YG1, R1

5 colors: + YG5


6 colors: +V6

7 colors: +V3

8 colors: +R3


+8 colors: Come to FACS Core Facility (you are welcome, also for fewer colors ☺)

Compensation Beads Positivity for each Fluorochrome. Analyzed on LSRFortessa 2015.05.20 BV586 analyzed 2016.02.08



Lymphocytes: Gating for calculating Stain Index

Analyzed on LSR Foreessa 2015.05.20 BV786 analyzed 2016.02.08




Stain Index for 25 different CD8 fluorochrome conjugates

| Fortessa, 200515, CD8 | Lasere | Filter | Fluorochrom | Stain Index | Median pos | Median neg | rSD neg |
|------------------------------|----------|--------|-----------------|-------------|------------|------------|---------|
| 200515_A1_A01_001.fcs | VIOLET | 442/46 | V450 | 46 | 5404 | 79,6 | 57,3 |
| 200515_A2_A02_002.fcs | VIOLET | 442/46 | BV421 | 370 | 80682 | 61,6 | 109 |
| 200515_A3_A03_003.fcs | VIOLET | 525/50 | V500 | 29 | 4682 | 186 | 77,9 |
| 200515_A4_A04_004.fcs | VIOLET | 525/51 | BV510 | 141 | 26259 | 166 | 92,8 |
| 200515_A5_A05_005.fcs | VIOLET | 610/20 | eVolve605 | 46 | 18458 | 277 | 198 |
| 200515_A6_A06_006.fcs | VIOLET | 610/20 | BV605 | 145 | 110959 | 310 | 382 |
| 200515_A7_A07_007.fcs | VIOLET | 660/20 | BV650 | 166 | 108268 | 198 | 326 |
| 200515_A8_A08_008.fcs | VIOLET | 711/25 | BV711 | 143 | 81675 | 162 | 286 |
| MNC_E3_E03_002.fcs | VIOLET | 785/60 | BV786 | 247 | 101078 | 123 | 204 |
| | | | | | | | |
| 200515_A10_A10_010.fcs | BLÅ | 530/30 | FITC | 53 | 9795 | 170 | 91,4 |
| 200515_A11_A11_011.fcs | BLÅ | 530/30 | AF488 | 82 | 16344 | 177 | 98,2 |
| 200515_A12_A12_012.fcs | BLÅ | 530/30 | BB515 | 150 | 48255 | 174 | 160 |
| 200515_B1_B01_013.fcs | BLÅ | 710/50 | PerCP-Cy5.5 | 121 | 18986 | 123 | 77,8 |
| 200515_B2_B02_014.fcs | BLÅ | 710/50 | PerCP-eFlour710 | 427 | 100831 | 82,2 | 118 |
| | | | | | | | |
| 200515_B3_B03_015.fcs | GUL-GRØN | 586/15 | PE | 537 | 102577 | 30,8 | 95,4 |
| 200515_B4_B04_016.fcs | GUL-GRØN | 615/20 | PE-eFlour610 | 342 | 86829 | 7,7 | 127 |
| 200515_B5_B05_017.fcs | GUL-GRØN | 615/20 | PE-CF594 | 327 | 113161 | 8,98 | 173 |
| 200515_B6_B06_018.fcs | GUL-GRØN | 660/20 | PE-Cy5 | 195 | 69002 | 47,5 | 177 |
| 200515_B7_B07_019.fcs | GUL-GRØN | 780/60 | PE-Cy7 | 266 | 94835 | 35,9 | 178 |
| | | | | | | | |
| 200515_B8_B08_020.fcs | RØD | 670/30 | APC | 169 | 38075 | 123 | 112 |
| 200515_B9_B09_021.fcs | RØD | 670/30 | AF647 | 208 | 63821 | 123 | 153 |
| 200515_B10_B10_022.fcs | RØD | 730/45 | AF700 | 66 | 22716 | 62,9 | 171 |
| 200515_B11_B11_023.fcs | RØD | 780/60 | APC-eFlour780 | 90 | 61082 | 87,4 | 337 |
| 200515_B12_B12_024.fcs | RØD | 780/60 | APC-H7 | 64 | 58606 | 146 | 458 |
| 200515_C1_C01_025.fcs | RØD | 780/60 | APC-Cy7 | 74 | 76083 | 130 | 514 |


Lasers and filters at the LSRFortessa

May 2016, Facs Core Facility, Aarhus University


405 nm laser Violet



488 nm laser Blue



561 nm laser Yellow-Green



640 nm laser Red

